

品名: Anti-Beta-actin Antibody (Mouse)

货号: AYC32-100 规格: 100ul

➤产品信息:

This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.

抗体来源: 小鼠

克隆类型: 单克隆抗体

来源: Recombinant protein of human β -actin 免疫

分子量: 41KD

适应种属: Mouse, Rat, Human, Fish (predicted_ Horse, Chicken, Guinea pig, Hamster, Cat, Dog, Pig, Xenopus laevis, Monkey, Zebrafish, Baboon, Xenopus tropicalis)

纯化方式: rProtein G 亲和纯化

产品纯度: >95%, SDS-PAGE 检测

储存缓冲液: 1xPBS (pH7.4), 1%BSA, 0.02%叠氮钠, 50%甘油

储存条件: 蓝冰运输, -20° C 可保存不少于 2 年, 避免反复冻融

➤产品应用:


WB=1/5000-1/20000

➤保存条件:

蓝冰运输, -80° C 存储>4 年, -20° C 可保存>2 年, 避免反复冻融

➤数据展示:

WB 数据:



Western blot analysis in various tissues and cell lines.